


ONLINE TEMPERATURE MONITORING




ADVANTAGES



- easy installation by user
- no investment costs, just rental costs
- automatic documentation and data backup
- predefined PDF reports by e-mail daily, weekly or monthly
- 6 years history available
- free replacement in case of defect
- other types of sensors available
- service and support contract included


ICUtechWebLog



The temperature monitoring is done via wireless sensor, which communicates with the ICU gateway. Measuring point basing on its user-defined limits. If a reading is outside of limits, an alarm is sent by email and/or SMS to a user. If the user does not handle the alarm, it is sent again an email and/or SMS after 24 hours.

Web interface

The web interface and the free app shows alarm status and temperature charts and the user can confirm the alarms and generate reports.



ICU tech WebLog is a product, which makes your work easier and gives you the quality you need.